

Gharieni Lina Select Podo Static

With its one-pillar technology and solid floor installation, the Gharieni Lina Select Static looks particularly elegant. At the same time it gives you practical legroom with maximum stability. In addition, this model has the same features as the "Wood" and "Alu."

Features:

- Bolted to the floor for maximum stability and legroom
- Four-actuator
- Pull-out and removable head and foot sections
- 5-actuator version: height, seat inclination, back and leg sections electrically adjustable using hand switch
- 2-actuator version: height and seat inclination electrically adjustable by foot switch, back and leg sections by gas spring
- Equipped with memory option, up to eight reclining positions can be stored with 5-actuator model
- Fold-down armrests upholstered on both sides
- Vertically adjustable foot section
- Upholstery with a memory foam layer
- Upholstery sewn into one piece
- "Airback" lordosis support system integrated into the upholstery (see page 20)



Lina Select Podo Static
with standard upholstery

Lina Select Podo Static

2-actuator
5-actuator

from € 2,971.60
from € 3,414.35



Lina Select Podo Static
with soft upholstery



- Height: approx. 51-93 cm
- Max. working height at the heel: approx. 125 cm
- Bed surface length: max. 202 cm
- Width without/with armrests: 62/80 cm
- (Recommended) lifting capacity: 200 kg
- Load capacity of the seat: 250 kg