

# Cryosurgical Systems: Fast, Precision Freeze and Thaw Control

**OUR ONE- AND TWO-TRIGGER SYSTEMS PROVIDE PRECISE CONTROL AND A LARGE SELECTION OF HIGHLY CONDUCTIVE, REUSABLE CRYOTIPS AND DISPOSABLE PLASTIC SHIELDS ARE AVAILABLE**



## LL100™ System (N<sub>2</sub>O)

Two-trigger Design and Twice the Freeze Power of Other Systems

Wallach's exclusive two-trigger cryosurgical system delivers fast and precise freeze/thaw control and allows cryotips to be changed at any time without venting.

- More effective freeze/thaw in less time
- Freeze and thaw only when triggered
- Lightweight, well-balanced, and maneuverable
- Use with Wallach's large selection of cryotips and disposable plastic shields

## LLCO<sub>2</sub>™ System

Single-trigger Design is Optimized for CO<sub>2</sub> Gas

The LLCO<sub>2</sub> single-trigger cryosurgical system maximizes temperature and minimizes the clogging and sputtering commonly associated with CO<sub>2</sub> gas.

- Squeeze trigger to freeze; release to thaw
- Use with medical-grade CO<sub>2</sub> gas for optimal performance
- Use with Wallach's standard cryotips with disposable plastic shields

Both the LL100 and the LLCO<sub>2</sub> are available in three-tip package specials for additional savings. Choose three cryotips from our comprehensive offerings when ordering. Wallach recommends the T-1905, T-2514, and T-0500 for GYN procedures and the T-0524, T-0219, and T-1900 for general/family practice procedures. Call Wallach customer service for more information.


ITEM NUMBER	DESCRIPTION
900001	LL100 Cryosurgical System – N <sub>2</sub> O (Freezer Only) <sup>Δ†</sup>
900160	LLCO <sub>2</sub> Cryosurgical System – CO <sub>2</sub> (Freezer Only) <sup>Δ</sup>
901060	N <sub>2</sub> O Cylinder (unfilled)
901061	CO <sub>2</sub> Cylinder (unfilled)
901011	Cart for 20 lb. Cylinder



## WA1000B™ Cryo Console System

- Uses LL100™, two-trigger or LLCO<sub>2</sub>™ single-trigger freezers
- High-performance features and easy-to-use controls
- Temperature Indicator with integrated timer for precise documentation
- Cart carries a single 20 lb. cylinder
- Use with Wallach's large selection of cryotips and shields

ITEM NUMBER	DESCRIPTION
900506-1-N2O	WA1000B Cryosurgical Console – Blank N <sub>2</sub> O LL100 <sup>Δ†</sup>
900506-1-CO2	WA1000B Cryosurgical Console – Blank CO <sub>2</sub> LL100 <sup>†</sup>
900506-2-N2O	WA1000B Cryosurgical Console – Digital Temp. (N <sub>2</sub> O) Indicator, Timer, Base, and LL100 <sup>Δ</sup>
900506-2-CO2	WA1000B Cryosurgical Console – Digital Temp. (CO <sub>2</sub> ) Indicator, Timer, Base, and LLCO <sub>2</sub> <sup>Δ</sup>



## Cryosurgical Tips for Use with LL100™ and LLCO<sub>2</sub>™

Highly conductive, thermal, reusable, sterilizable cryotips with disposable plastic shields.

900300AA (T-0200) <sup>Δ†</sup>	900301A (T-0500) <sup>Δ†</sup>	900302A (T-0700) <sup>Δ†</sup>	900303AA (T-0900) <sup>Δ†</sup>	900305AA (T-1300) <sup>Δ†</sup>	900306AA (T-1500) <sup>Δ†</sup>
900200AA (T-0219) <sup>Δ†</sup>	900201AA (T-0219) <sup>Δ†</sup>	900202AA (T-0524) <sup>Δ†</sup>	900203AA (T-0519) <sup>Δ†</sup>	900204AA (T-0832) <sup>Δ†</sup>	900205AA (T-0819) <sup>Δ†</sup>
900206AA (T-0826) <sup>Δ†</sup>	900207AA (T-0823) <sup>Δ†</sup>	900208AA (T-1900) <sup>Δ†</sup>	900209AA (T-1905) <sup>Δ†</sup>	900210AA (T-1920) <sup>Δ†</sup>	900211AA (T-1910) <sup>Δ†</sup>
900212AA (T-1920) <sup>Δ†</sup>	900213AA (T-2500) <sup>Δ†</sup>	900214AA (T-2507) <sup>Δ†</sup>	900215AA (T-2514) <sup>Δ†</sup>	900217AA (T-0832) <sup>Δ†</sup>	900221AA <sup>Δ†</sup>

Disposable Plastic Shields for Cryotips (Bag of 100)

