

A COMPACT SOLUTION THAT ERADICATES BACTERIA, VIRUSES, & FUNGI

ETO STERILIZER FULLY AUTOMATIC



PLT-103(I)



- Versatile solution for sterilizing medical devices, from single-use items to surgical instruments.
- It is used in Surgical threads, Surgical trays Procedure kits Synthetic gowns Fine surgical instruments
- ETO can kill all known viruses, bacteria and fungi, including bacterial spores.
- Ethylene oxide treatment is generally carried out between 30°C and 60°C with relative humidity above 30%.
- Display of set values, temperature, pressure & humidity on HMI screen.
- Manufactured under GMP norms, ensuring compliance with industry standards



MODEL	ETOS-1	ETOS-2	ETOS-3	ETOS-4
Working Volume(Ltrs.)	55 Ltrs.	81 Ltrs.	110 Ltrs.	180 Ltrs.
Working Capacity(Cu. Ft.)	2 Cu. Ft.	2.8 Cu. Ft.	4 Cu. Ft.	6.7 Cu Ft.
Chamber Size	12 x 12 x 24	12 x 12 x 36	12 x 12 x 48	18 x 18 x 36
Sterilant Gas	40 grams cartridge	100 grams cartridge	100 grams cartridge	100 grams cartridge



P.L.TANDON & CO.

Office: 205-B Apra Plaza Building Road No:44,
Pitampura, New Delhi - 110034

Factory: I-59, South City -2, Sector - 40,
Gurugram, Haryana - 122004

Contact: +91-9810347149, +91-9999054454

info@tancolabproducts.com www.tancolabproducts.com

