



LIFETIME SA01PT



Product Page

Is a suction unit used for body fluids, oral, nasal and tracheal aspirations in adults and children.



Compact Design

Thanks to its ergonomic design, it provides ease of use and storage.



Adjustment Knob

You can reach the desired pressure in the shortest time with the sensitive adjusting knob.



Autoclavable Jar

Thanks to its autoclavable polycarbonate jars at 121 degrees, it provides hygienic use for many times.



Oil Free Piston Pump

Its special compressor does not require maintenance and allows use for many years.

Technical Specifications

Compressor	Oilless Piston Pump
Voltage	230 VAC \pm %10 , 50 Hz
Max. Free Air Flow	12 lt / min \pm %10
Max. Vacuum	-550 mmHg \pm %5
Power Consumption	0,60 A
Vacuum Setting	Yes
Collection Jar	1 lt
Float System	Yes
Hydrophobic Filter	Yes
Standard	ISO 10079 – 1 / CE
Device Dimensions	20x15x25 cm
Weight of Device	2.5 kg
Box Dimensions	34x16x28 cm
Gross Weight (Device)	2.8 kg
Stock Code	PT01B000
Warranty Period	2 years



20-30
LPM

LIFETIME SA03PB/SA04PB



Product Page

Ideal for ambulances, emergency services, examination rooms, inpatient services, outpatient clinics and home care.



Analog Display

Its analog screen provides the most accurate way to monitor the adjusted pressure with a sensitive vacuum meter.



Adjustment Knob

You can reach the desired pressure in the shortest time with the sensitive adjustment knob.



Li-ion Battery

Provides up to 60 minutes of use on a single full charge. No need for battery replacement for a long time.



Wall Bracket

Thanks to its wall bracket that comes as a standard in clinics and ambulances, stable use is ensured.



LED Charge Indicator

Thanks to the LED charge indicator, the battery charge rate can be monitored continuously in 5 stages.



Accessory Box

The accessory box provides easy transportation and use of all consumables with the device, especially for mobile uses.

Technical Specifications

Product Model	SA03PB	SA04PB
Flow Rate	20 LPM \pm %10	30 LPM \pm %10
Noise Level	<50 dB	<60 dB
Weight	4 kg	4.5 kg
Battery	90 min	60 min
Max. Vacuum	-650 \pm 5 mmHg	
Power Consumption	30W	60W
Compressor Type	Oilless Piston Pump	
Battery Type	Li-ion	
Voltage	12 VDC	
Jar Capacity	1 lt	
Box Dimensions	40x20x33 cm	
Gross Weight	5 kg	6 kg
Warranty Period	2 years	
Stock Code	PB03B000	PB04B000