

# Ameda® Pearl™ Hospital Grade Breast Pump



Trusted by Hospitals, Loved by Moms.

## Powerful Performance in the Hospital and Home.

Combines powerful performance inside with a protective shell on the outside. Pearl is designed to build milk supply amongst pump dependent, late pre-term and term moms. Pearl's modern design combines the durability needed for the multi-user/hospital environment with lifestyle features moms expect. Pearl will serve the needs of hospitals and is the ideal multi-user product for milk banks and rental stations. Pearl is perfectly positioned to support continuity of care in the Hospital and Home.

### Ameda Pearl – Proven Hospital Grade Pumping

- **Proven:** Over 90% of moms surveyed were highly satisfied with Pearl
- **Powerful:** Hospital Grade / Multi-User Performance with two modes of pumping
- **Portable:** Lightweight with rechargeable battery and integrated handle
- **Smart:** Vacusense™ Technology ensures consistent vacuum performance in varying altitudes or changing atmospheric conditions

- **Comfortable:** Fully customizable expression mode with 72 possible combinations of speed and suction settings
- **Durable:** Protective housing designed for easy cleaning between users; 3-year limited warranty
- **Whisper Quiet:** Discreet pumping anywhere
- **Safe:** Compatible with Ameda's Universal HygieniKit™, the only FDA-cleared closed system pump kit proven to protect against contaminants like bacteria, viruses and mold growth



### Ameda Pearl Technical Specifications

Suction Range	Stimulation: 50 - 150mmHg; 6 user-selectable settings Expression: 30 - 250mmHg; 12 user-selectable settings
Cycles per Minute (CPM)	Stimulation: 80 - 120 CPM; 3 preprogrammed settings Expression: 30 - 48 CPM; 6 user-selectable settings
Power Supply	Medical AC Adapter – Output: 15VDC, 2A, 30W Input: 100-240VAC, 50/60Hz, 1A 2600mAh Lithium-ion rechargeable battery; 3 hours per charge
Unit Weight	1.6kg
Unit Dimensions	20 x 18 x 15cm (L x W x H)
Preprogrammed pumping	Mode auto switches from stimulation to expression at 2 mins Auto Power off at 60 mins Vacusense Technology
Trolley	Optional

Authorised Distributors :

**JOLLY TECHNOCRATS**

An ISO 9001:2015 Company

Complete Range of : OT Equipment, Skin, Surgery & Gynae Products

304, Mansarovar Building, 90, Nehru Place, New Delhi-110019, (INDIA)

Phone : + 91-11-46066071, Telefax : + 91-11-46066072,

e-mail : jollytechnocrats@gmail.com, info@jollytechnocrats.in, myamedain@gmail.com

website : www.jollytechnocrats.in, Mobile : +91-9810520055, +91-9650595508



Ameda Pearl Trolley  
201H01

