



MOM INSPIRED. HOSPITAL TRUSTED.

Personal Breast Pumps

Personal Breast Pumps

The simple and effective choice for mothers returning to work or who want a quick, easy and effective breast pumping experience.

Lactaline™ Double Electric Breast Pump

- Separate suction and speed dials allow for CustomControl™
Dual adjustability lets the mother select the best levels for her milk flow and comfort, creating a multi-phase pumping experience
- Airlock Protection™ prevents airflow between expressed milk and pump tubing while pumping
- Convenient power source options:
AC adapter, AA batteries or optional car adapter

Includes:

- Lactaline Electric Breast Pump
- Double Hygiene Pumpset
- AC adapter



Pumpset
without BPA

One-Hand Breast Pump

with Flexishield™ Areola Stimulator

- Ergonomic and compact manual pump earned the High Design-Innovation Quality Award from Design Zentrum Nordrhein-Westfalen Germany
- Pump directly into any standard baby bottle or freezer bag
- Allows mothers to achieve multi-phase pumping by varying handle squeezes for more milk flow
- Pump includes the Flexishield Areola Stimulator to help stimulate the let-down reflex



No BPA

Personal Breast Pumps

Power Adapters



AC Adapter: Europe




AC Adapter: UK



AC Adapter: Australia



12V Car Adapter




Lactaline Double Electric Breast Pump




One-Hand Breast Pump with Flexishield


Model Number	17080EA	17080EB	17080EC	17080UK	17066I
Power Supply	EU	EU	EU	UK	n/a
Minimum Order Qty.	1 each	1 each	1 each	1 each	1 each
Contents per Min. Order	1 piece	1 piece	1 piece	1 piece	1 piece
Product Dimensions (cm)	14.0L x 22.9W x 22.2H	14.0L x 22.9W x 22.2H	14.0L x 22.9W x 22.2H	14.0L x 22.9W x 22.2H	1.7L x 2.5W x 3.5H
Product Weight (kg)	1.8	1.8	1.8	1.8	0.3
Case Dimensions (cm)	15.2L x 25.4W x 25.4H	15.2L x 25.4W x 25.4H	15.2L x 25.4W x 25.4H	15.2L x 25.4W x 25.4H	39.4L x 52.7W x 21.6H
Case Weight (kg)	1.8	1.8	1.8	1.8	3.8
Qty. per Carton	n/a	n/a	n/a	n/a	12 each
Packaging Type	Retail Box DE - FR - IT - NL	Retail Box SE - DK - NO - FI	Retail Box EN - ES - PT - GR	Retail Box EN	Retail Box




AC Adapter EU
230V, 50 Hz/9v, 1000mA



AC Adapter UK
230V, 50 Hz/9v, 1000mA



AC Adapter AU
230V, 50 Hz/9v, 1000mA



Car Adapter
12V



Lactaline
Base Unit Only

Model Number	402706	402705	402707	17079I	2450208I
Minimum Order Qty.	1 each	1 each	1 each	1 each	1 each
Contents per Min. Order	1 piece	1 piece	1 piece	1 piece	1 piece
Product Dimensions (cm)	5.0L x 5.0W x 7.7H	5.0L x 5.0W x 5.0H	5.0L x 5.0W x 7.7H	2.5L x 2.5W x 7.6H	14.0L x 22.9W x 22.2H
Product Weight (kg)	0.46	0.46	0.46	0.11	0.64
Case Dimensions (cm)	45.7L x 30.5W x 61.0H	45.7L x 30.5W x 61.0H	40.0L x 32.0W x 18.0H	27.9L x 30.5W x 45.7H	56.0L x 39.8W x 20.8H
Case Weight (kg)	16.3	18.1	19.0	17.0	8.7
Qty. per Carton	36 each	40 each	40 each	150 each	12 each
Packaging Type	Bulk Box	Bulk Box	Bulk Box	Retail Box	Bulk Box