



MOM INSPIRED.HOSPITAL TRUSTED.

Ameda Una single electric breast pump, designed for daily use, is the ideal breast pump.

Single electric dual Phase Breast Pump. Features of dual phase Control like a hospital pump. Long term use helps in improving expressions as well as stimulate milk.

Includes:

- ❖ Single Electric Breast Pump, dual Phase
- ❖ Single Eco Pump set
- ❖ AC Adapter

The Ameda Difference

Purity.

Ameda has the world's only breast pump kit with a proven protective barrier. Proven Airlock Protection™ keeps tubing dry by preventing air flow between the expressed milk and pump tubing while pumping. No need to clean the narrow tubing and milk stays pure.

Comfort.

The Ameda CustomFit Flange System™ makes it easy to find a breast flange fit by offering seven breast flange sizes that adapt to all Ameda breast pumps and kits (available separately).

Ease.

Ameda breast pumps feature fully independent CustomControl™ settings for speed and suction. Dual adjustability lets the mother select her best levels for milk flow and comfort, so that she can create a multi-phase pumping experience.



Single Eco pumpset kit

Includes:

- ❖ 1 Breast Shell
- ❖ 1+1 Valves (spare)
- ❖ 1 diaphragm
- ❖ 1 Bottle, 1 tubing
- ❖ 1 cap, 1 Bottlecap
- ❖ 1 Adapter



For Sale & also available on Rental

Authorised Distributors :

JOLLY TECHNOCRATS

Complete Range of : OT Equipment, Skin, Surgery & Gynae Products

304, Mansarovar Building, 90, Nehru Place, New Delhi-110019, (INDIA)

Phone : + 91-11-46066071, Telefax : + 91-11-46066072,

e-mail : jollytechnocrats@gmail.com, info@jollytechnocrats.com

website : www.jollytechnocrats.com, Mobile : +91-9810520055, +91- 9650595508



Personal Single Electric Breast Pump