


CustomFit® Flange System

The Ameda CustomFit Flange System makes it easy to find a good flange fit by offering seven sizes that adapt to all Ameda Breast Pumps and Kits. Research has found that nearly one third of mothers who pump report greater comfort and more milk with a larger-than-standard breast flange size.^{1,2} Ameda CustomFit® Flanges have been designed to work with any Ameda HygieniKit® Milk Collection Kit. Simply transfer the valves and diaphragms to the new, better-fitting breast flange size and you are ready to pump.



Flexishield™ Areola Stimulator




Reducing Insert



- XS** Flexishield™ Areola Stimulator | 21.0mm
(Fits Inside Standard 25.0mm Breast Flange)
- S** Reducing Insert | 22.5mm
(Fits Inside Standard 25.0mm Breast Flange)

Standard Breast Flange



STANDARD

Standard Breast Flange | 25.0mm
(Included with all Ameda HygieniKit Systems and Breast Pumps)

L • Custom Breast Flange



M • Custom Breast Flange



- M** Custom Breast Flange | 28.5mm
- L** Custom Breast Flange | 30.5mm

XXL • Custom Breast Flange



XL • Custom Breast Flange Insert



- XL** Custom Breast Flange Insert | 32.5mm
(Fits Inside XX-Large 36.0mm Breast Flange)
- XXL** Custom Breast Flange | 36mm

¹ Jones, E. et al. A randomized controlled trial to compare methods of milk expression after preterm delivery. Arch Dis Child Fetal Neonatal Ed 2001; R91-F95.

² Wilson-Clay, B. and K. Hoover. The Breastfeeding Atlas, 4th ed. Austin, TX: LactNews Press, 2008.

CustomFit Flanges

Provides mothers with larger flanges and special reducing inserts to create an “in-between” size.

- ~ (2) Breast flanges
- ~ (2) Inserts
- ~ Without BPA or DEHP



Item Number	Description	Sale Unit	UPC
17148IP*	L/M (30.5mm flange/ 28.5mm insert)	1 case (6 pair)	810725020975
17358I*	XXL/XL (36.0mm flange/ 32.5mm insert)	1 case (6 pair)	810725021149

*Must be ordered in increments of 10

CustomFit Flange System Spare Parts

Item Number	Description	Sale Unit	UPC
10011*	22.5mm Reducing Insert	1 each	810725023129
10012	25.0mm Standard Breast Flange	1 each	810725023143
24110725	25.0mm Breast Flange – Bulk	1 each	NA
10013*	28.5mm Breast Flange	1 each	810725023150
10014*	30.5mm Breast Flange	1 each	810725023167

*Must be ordered in increments of 10

Flexishield™ Areola Stimulator

Gentle massaging action, softens the feel of the breast flange, reduces standard 25.0mm flange size to 21.0mm for maximum effectiveness for smaller nipples.

- ~ Without BPA or DEHP




Item Number	Sale Unit	UPC
17212I	1 pair	810725021002
10010*	1 each	810725023112

*Must be ordered in increments of 10


Did you know?

Incorrect flange fit can affect breast pump suction! Using the right size breast flange can help make pumping more comfortable and efficient.




Good Fit

During pumping, your nipple moves freely in the breast flange tunnel. You see space around the nipple. Not much areola is drawn into the tunnel with the nipple.



Too Small

During pumping, some or all of your nipple rubs against the sides of the breast flange tunnel.



Too Large

During pumping, more areola is drawn into the breast flange with your nipple. Your areola may rub against the side of the breast flange tunnel.