

Universal HygieniKit Options



PRODUCT	SKU	UPC	Qty*	Unit Dimensions	GTIN
STERILE					
Universal Single HygieniKit w/25mm flange 2 valves; 1 diaphragm; 1 tubing; 1, 4oz/118ml bottle; 1 tubing adapter & cap; 1, 25mm flange; 1 locking disc & ring	200H02	810725024102	18	36 x 20 x 10cm (H x W x D)	10810725024109
Universal Single HygieniKit w/28.5mm flange 2 valves; 1 diaphragm; 1 tubing; 1, 4oz/118ml bottle; 1 tubing adapter & cap; 1, 28.5mm flange; 1 locking disc & ring	200H03	810725024096	18	36 x 20 x 10cm (H x W x D)	10810725024093
Universal Single HygieniKit w/30.5mm flange 2 valves; 1 diaphragm; 1 tubing; 1, 4oz/118ml bottle; 1 tubing adapter & cap; 1, 30.5mm flange; 1 locking disc & ring	200H04	810725024140	18	36 x 20 x 10cm (H x W x D)	20810725024144
Universal DUAL HygieniKit w/25mm flange 4 valves; 2 diaphragms; 2 tubing; 2, 4oz/118ml bottles; 1 tubing adapter; 2 tubing adapter caps; 2, 25mm flanges; 2 locking discs & rings	200H05	810725024119	16	38 x 25 x 10cm (H x W x D)	10810725024116
NON-STERILE					
Universal Single HygieniKit w/25mm flange 2 valves; 1 diaphragm; 1 tubing; 1, 4oz/118ml bottle; 1 tubing adapter & cap; 1, 25mm flange; 1 locking disc & ring	200U02	810725025949	18	36 x 20 x 10cm (H x W x D)	10810725025946
Universal Single HygieniKit w/28.5mm flange 2 valves; 1 diaphragm; 1 tubing; 1, 4oz/118ml bottle; 1 tubing adapter & cap; 1, 28.5mm flange; 1 locking disc & ring	200U03	810725025956	18	36 x 20 x 10cm (H x W x D)	10810725025953
Universal Single HygieniKit w/30.5mm flange 2 valves; 1 diaphragm; 1 tubing; 1, 4oz/118ml bottle; 1 tubing adapter & cap; 1, 30.5mm flange; 1 locking disc & ring	200U04	810725026069	18	36 x 20 x 10cm (H x W x D)	10810725026066
Universal DUAL HygieniKit w/25mm flange 4 valves; 2 diaphragms; 2 tubing; 2, 4oz/118ml bottles; 1 tubing adapter; 2 tubing adapter caps; 2, 25mm flanges; 2 locking discs & rings	200U05	810725026076	16	38 x 25 x 10cm (H x W x D)	10810725026073

Single Eco Pumpset

Sterile

- ~ (1) Tubing adapter
- ~ (1) Silicone tubing
- ~ (1) Adapter cap
- ~ (1) Silicone diaphragm
- ~ (1) 25.0mm breast flange
- ~ (1) Valve
- ~ (1) 4 oz. bottle with solid sealing lid
- ~ Without BPA or DEHP

Item Number	Sale Unit	UPC
171051*	1 case/24 each	810725022412

*Must be ordered per case



Single HygieniKit Pumpset

Non-Sterile

- ~ (1) Tubing adapter
- ~ (1) Silicone tubing
- ~ (1) Adapter cap
- ~ (1) Silicone diaphragm
- ~ (1) 25.0mm breast flange
- ~ (2) Valves
- ~ (1) 4 oz. bottle with lock-tight sealing lid
- ~ Without BPA or DEHP

While
Supplies
Last

Item Number	Sale Unit	UPC
171111P*	1 each	810725021156

*Must be ordered in increments of 20

