



Trusted by Hospitals, Loved by Moms.

AMEDA GLO™ Wearable Breast Pump



What's included



24mm FLANGES (2)



DIAPHRAGMS (2)



VALVES (4)



WEARABLE PUMPS (2)



17, 19, & 21mm
FLANGE INSERTS (2 EA)







USB-C CHARGING CABLES (2)



BRA ADJUSTMENT STRAPS (2)

Meet Ameda GLO – the wearable pump for modern moms

-  Hospital grade performance
-  Milk optimization technology
-  4 pumping modes & 15 levels of suction
-  Pumping level memory



Travel case sold separately

SKU #	GLO PUMP & ACCESSORIES
120G01	Your Ameda Glo Includes: (2) Glo Wearable pumps; (2) 6 oz 177mL collection cups; (2) USB charging cables; (2) 24mm flanges; (2) 17mm, (2) 19mm & (2) 21mm soft flange inserts; IFU; (2) bra extenders; (4) valves
120G02	Ameda Universal Breast Pump Travel Case (wearable pump sold separately)
120G03	Ameda Glo 24 mm flange, non-sterile, 1 ea.
120G04	Ameda Glo 27 mm flange, non-sterile, 1 ea.
120G05	Ameda Glo Diaphragm, 1 ea.
120G06	Ameda Glo Valves, 2 ea.
120G07	Ameda Glo USB Charging Cable, 1 ea.
120G08	Ameda Glo Bra Adjustment Buckles, 2 ea.
120G15	Ameda Glo Collection Cup, 1 ea.

SKU #	AMEDA UNIVERSAL SOFT FLANGE INSERTS
120G09	Ameda Universal Soft Flange Inserts 11mm non-sterile, 2 ea.
120G10	Ameda Universal Soft Flange Inserts 13mm non-sterile, 2 ea.
120G11	Ameda Universal Soft Flange Inserts 15mm non-sterile, 2 ea.
120G12	Ameda Universal Soft Flange Inserts 17mm non-sterile, 2 ea.
120G13	Ameda Universal Soft Flange Inserts 19mm non-sterile, 2 ea.
120G14	Ameda Universal Soft Flange Inserts 21mm non-sterile, 2 ea.
120G16	Ameda Soft Flange Inserts Multipack 11mm, 13mm, 15mm, 17mm, 19mm, 21mm non sterile, 2 ea.



SPECIFICATIONS	
Dimensions & Weight	4.88" x 4.37" x 2.68" 124.9mm x 111mm x 68.1mm; Weight: 8.6 oz 243.8g (each)
Collection Cup Capacity	6 oz 177mL
Modes	Stimulation, Expression, Massage, Auto



Trusted by Hospitals, Loved by Moms.

Ameda.com | ParentCare@ameda.com | 1.866.992.6332

