



LOOKS[®] SURGICAL TELESCOPES

Manufacturer : **XENOSYS CO., LTD.**

2, Venture-ro, 100 beon-gil, Yeonsu-gu, Incheon 406 840, Republic of Korea

Tel. +82 32 875 9811 Fax. +82 32 875 9813

E-mail : sales@xenosys.co.kr / hadlee@hanmir.com URL : <http://www.xenosys.co.kr>

FEATURES

- Larger field of view & superior depth of focus
- Lighter & Smaller Telescopes
- Clear & Bright vision over the full field
- No Distortion image
- No Fatigue on eye
- IP Frames
 - Super light weight & comfortable aircraft aluminum front
 - Hypoallergenic titanium temple
 - Ion Plating Color

• IPHL Frame



• IPHX Frame



SPECIFICATIONS

| Model | Mag. | F.O.V. @40cm | Weight* |
|-----------|---|----------------|---------|
| Looks2500 | 2.5x | 160 mm @ 40 cm | 33.0g |
| Looks2800 | 2.8x | 130 mm @ 40 cm | 31.8g |
| Looks3200 | 3.2x | 86 mm @ 40 cm | 49.1g |
| Looks3500 | 3.5x | 98 mm @ 40 cm | 69.1g |
| Looks3700 | 3.7x | 85 mm @ 50 cm | 71.0g |
| Looks5000 | 5.0x | 62 mm @ 40 cm | 77.4g |
| Looks7500 | 7.5x | 41 mm @ 40 cm | 83.6g |
| Flip-up | Magnification : 2.5x Working Distance choice among: 30 cm, 36 cm, 42 cm, 48 cm | | |

* Weight based on TTL loupe with plano lens with IPAL frame

Flip-up 2.5x



Looks 2500



Looks 2800



Looks 3200



Looks 3500



Looks 3700



Looks 5000



Looks 7500



FDA CE CFDA ISO13485:2016 ISO14001:2015



Authorised Distributors :

JOLLY TECHNOCRATS

An ISO 9001:2015 Company

Complete Range of : OT Equipment, Skin, Surgery & Gynae Products

304, Mansarovar Building, 90, Nehru Place, New Delhi-110019, (INDIA)

Phone : + 91-11-46066071, Telefax : + 91-11-46066072, e-mail : jollytechnocrats@gmail.com, info@jollytechnocrats.in

website : www.jollytechnocrats.in, Mobile : +91-9810025581, +91-9650595509