



The Whole brand-new LED Headlight System

▪ L2SX5 (Single or Dual LED choices)

- ✓ Maximum Brightness 110,000 lx. (LLMX566D model, Battery powered @40cm)
- ✓ Single LED X566 have 61,000 lx. (Distance @ 40 cm away)
- ✓ Magnetic extension cable: Easy to attach & detach, comfortable to store
- ✓ Comes with light-weight handy carrying case ▪ Battery running time maximum 24 hours (at the lowest level)

New Launch

Additional battery (ONCLUDED)



Remote Controller



Carrying Case



L2SX5-S (Single LED)				
Step	Brightness (lx) @WD400mm Illumination (±10%)			Operating Time
	LLMX582	LLMX574	LLMX566	
1 st	12,000	14,000	18,000	24H
2 nd	22,000	25,000	33,000	13H
3 rd	30,000	40,000	44,000	10H
4 th	41,000	47,000	61,000	6h30m
L2SX5 -D(Double LED)				
Step	Brightness (lx) @WD400mm Illumination (±10%)		Operating Time	
	LLMX566D			
1 st	36,000		12h50m	
2 nd	64,000		7h30m	
3 rd	87,000		5h00m	
4 th	110,000		3h10m	

FEATURES

- Super bright Dual LED headlight even its battery powered with ONE control unit
- Compact & light weight control unit, LED module
- Improved & compact size battery pack & Secure assembling system
- LED heat conduction effect: Multi link & attachment is made with aluminum
- LED module weight 12.3 g (LLMX582)
- More soft & durable LED wires



www.xenosys.co.kr

Manufacturer: XENOSYS CO., LTD.
 2, Venture-ro, 100 beon-gil, Yeosu-gu, Incheon 406-840, Republic of Korea
 Tel: +82-32-875-9811 Fax: +82-32-875-9813
 E-mail: sales@xenosys.co.kr / hadlee@hanmir.com URL: http://www.xenosys.co.kr