



**IMAGING EXCELLENCE  
OF 4K RESOLUTION**



- ◆ Resolution : 3840 x 2160
- ◆ Inputs Signals : 2x HDMI, 2x DP
- ◆ Output Signal : DP

- ◆ 27" MHSD-2740
- ◆ 32" MHSD-3240
- ◆ 43" MHSD-4340
- ◆ 55" MHSD-5540

- ◆ Enhanced Precision into Narrow & Wide Laparoscopic Surgical Procedures
- ◆ Touch Screen Optional
- ◆ PIP/ PBP (L/R, Top/Bottom) / 4P(4Windows)

- ◆ Clarity, Reliability for the Clear View of Video & Still Images from Endoscopes, Laparoscopic Cameras & Other Imaging Systems