



## Battery Type

DC 14,8 V, 26 Ah, 3,8 kg  
3 hours when fully charged(100%)



- Max 130,000 lx Intensity
- Intensity Control in 3 steps
- 3 Steps of color Temperature selectable (3,800K / 4,300K / 4,800K)
- Ra 96 of Highest Level of CRI
- Adjustable Focus (Min 15cm, Max 25cm)
- Excellent Dilution Effect with 42 LED Units
- No UV emission & Lower Level of heat emission
- Detachable & Autoclavable Hand Grip

