

Spider Arm

Effective Capable Easy



Robot arm for multi-function joints



Capable of overcoming high loads



High reliability and efficiency



Easy to use and maintain



Surgery Head Frame

Flexible Precise Support



Provides excellent flexibility




High support enables stable surgery



Precise movement can be implemented



Creates the optimal environment for nerve surgery



Size (mm)	950×300×80	Load capacity (kg)	250
Rotation angle (°)	360	Joints	6

Application	Neurosurgery	Rod diameter (mm)	12-20
Inserting rods distance (mm)	220-480	Package	Carton