

## Spine Surgery Frame

Comfortable Effective Convenient



Fast and convenient locker operation



Matches all brands of surgical tables



Improving the effectiveness of spinal surgery



Provide comfort to patients



Size (mm)	750×500×180-260	Material	Carbon fiber
Operation method	Bilateral manual crank operation	Color	Black

## Beach Chair

Flexible Suitable Convenient



Provides high flexibility



Convenient and quick surgical environment



Suitable for shoulder and beach surgery



Change the configuration according to the purpose of the surgery



Size (mm)	520×580	Back angle (°)	-15/80
Headrest adjustment Range (mm)	0-320	Required operating table Width (mm)	500-600